

AIROS® 6 SEQUENTIAL COMPRESSION DEVICE



Reliable, Intuitive 6-Chamber Device Ideal for the Medicare Population

The FDA-cleared AIROS® 6 Sequential Compression Device delivers compression therapy to help patients manage lymphatic and venous disorders. The device provides effective intermittent pneumatic compression (IPC) therapy for treatment of the lower and upper extremities.

The AIROS 6 provides gradient sequential compression that increases blood flow and encourages extracellular fluid clearance. The device provides cycles of compressed air at certain adjustable pressures and then inflates and deflates garments to move fluid.

KEY FEATURES

- **Regulatory Clearances:** FDA (510K Number K172770) & Medicare (HCPCS Code E0651); also reimbursable through commercial insurance
- **Usage Tracking:** Built-in timer records total usage time to measure patient compliance
- **Easy Operation:** Large, clear LCD screen and one-button operation for starting and stopping therapy
- **Physical Benefits:** Reduce pain and swelling, increase blood flow, move lymphatic fluid
- **Various Sizes of Comfortable Garments:** 8 full leg sizes, 2 half leg sizes, and 2 arm sizes available to comfortably fit patients (all utilize 6 individual chambers)
- **Lymphedema Preparation Treatment (LPT) Modes:** Optional Clearing Mode (mimics lymphatic drainage) or Focal Mode (fibrotic tissue relief) can be used before main Gradient Mode therapy



SPECIFICATIONS

Dimensions	10.2" W x 10.5" D x 4.3" H
Weight	6.4 lbs.
Primary Operating Mode	Gradient
Secondary Operating Modes (LPT)	Clearing or Focal
Pressure Range (Gradient Mode)	30-80mmHg (increments of 1mmHg)
Cycle Time (Inflation)	40 or 60 seconds
Deflation Time	12 seconds
Therapy Time	15, 30, 45 or 60 minutes
Warranty (Device)	3 years
Alarms (Error Messages)	3
Air Blocker (Unilateral Therapy)	Included (1)
Electrical	120 VAC, 60Hz, 100VA MAX