



**Questions?
Please Contact Us Immediately.**

Dedicated to Breast Cancer Patient Care

In the United States, approximately 266,000 new cases of breast cancer are diagnosed each year. Many of these patients will develop lymphedema. Our products, designed with substantial input from lymphedema therapists, provide upper extremity treatment that reduces pain and swelling. We strive to help breast cancer survivors live more comfortably.

Personalized Care at Home

Patients can utilize AIROS devices and garments inside the comfort of their own homes. A highly-trained field representative sets the device pressure and timing according to their script. Customer service teams are available to answer any calls as needed.

Medicare Coverage with No Trial Needed

Our advanced AIROS 6 device and associated 6-chamber Arm Plus garments are covered by Medicare. No trial is necessary to begin patient therapy.

Covered by Commercial Insurance

Most compression therapy devices are also covered by private insurance. Billing experts can work with clinicians' offices to complete necessary paperwork and clinical notes.



AIROS[®]
MEDICAL

**A New Solution
for Upper Extremity
Lymphedema Treatment**

www.airosmedical.com

Compression therapy devices utilized with easy-to-apply Arm Plus garments help patients experience relief from chronic pain and swelling related to lymphedema.

www.airosmedical.com

© AIROS[®] Medical, Inc. 2020

All Rights Reserved.



AIROS 6 Sequential Compression Device

How It Works

The digitally-controlled AIROS 6 device inflates air through hoses into garments worn on the arm or arm, shoulder, and chest areas. The garments inflate at a set pressure (according to physician's script) in a specific sequence, beginning with the first chamber near the hand. As the garments inflate, fluid is moved to reduce swelling and pain over time. This therapy is known as gradient sequential compression.

Device Features

- ✓ Advanced device for Medicare population (HCPCS E0651) or commercial insurance
- ✓ Large, clear LCD screen for high visibility
- ✓ Lightweight, quiet, and easily portable
- ✓ Countdown timer during therapy & compliance-driven Usage Data Timer
- ✓ Built-in vacuum aiding in garment application and removal
- ✓ 3 Operating Modes
 - ▶ Gradient Mode: Main gradient sequential operating mode.
 - ▶ Focal Mode: Secondary operating mode providing focused therapy on areas with fibrotic tissue build-up (3 chambers active)
 - ▶ Clearing Mode: Secondary operating mode mimicking Manual Lymphatic Drainage (MLD)

AIROS 8 Sequential Compression Device

How It Works

The digitally-controlled AIROS 8 utilizes the same type of gradient sequential compression therapy as the AIROS 6, but utilizes two additional air chambers. The AIROS 8 device also features an additional operating mode that allows for custom pressures in each garment, allowing the user to calibrate the pressure (if prescribed by the physician).

Device Features

- ✓ High-end device for private insurance, also Medicare reimbursable (HCPCS E0652)
- ✓ Large, clear LCD screen for high visibility
- ✓ Lightweight, quiet & easily portable
- ✓ Countdown timer during therapy & compliance-driven Usage Data Timer
- ✓ Built-in vacuum aiding in garment application and removal
- ✓ 4 Operating Modes
 - ▶ Gradient Mode: Main gradient sequential operating mode.
 - ▶ Pressure Mode: Main gradient sequential operating mode, but allows user to set pressure in individual chambers or "skip" pressure in a chamber
 - ▶ Focal Mode: Secondary operating mode providing focused therapy on areas with fibrotic tissue build-up (4 chambers active)
 - ▶ Clearing Mode: Secondary operating mode mimicking Manual Lymphatic Drainage (MLD)

AIROS Compression Garments

Carefully Designed Garments Sized for All Patients

AIROS compression garments are utilized with the AIROS 6 (6-chamber garments) and AIROS 8 (8-chamber garments). The garments can easily be plugged into the front of the device before therapy. Made from a durable material, the garments are designed to support patients of all sizes and therapy needs.

Device Features

- ✓ 30 total sizes available
 - ▶ 6-Chamber: 2 Arm, 4 Arm Plus (Arm/Chest/Shoulder), 8 Full Leg
 - ▶ 8-Chamber: 2 Arm, 4 Arm Plus (Arm/Chest/Shoulder), 8 Full Leg, 2 Half Leg
- ✓ Arm Plus garments available in two colors and packaged with 3 garment extension pieces for proper fitting
- ✓ Support padding in the garment and O-Rings on the zippers help make garment application and removal easy/comfortable, easy-to clean material ideal for in-home use